Hot-air sterilizers have been equipped with a couple of additional functions that protect samples. They can sterilize at temperatures of up to 250°C



Sterilizer SRW 240 Smart IG



All thermostatic equipment manufactured by POL-EKO-APARATURA can be provided with Calibration Certificate issued by accredited Measurement Laboratory. Detailed information on accreditation is available on website: www.pol-eko.eu.

Heating equipment Sterilizers SR



The highest level of customer service is our priority. Our goal is to be a consulting company. We do our best to ensure that our sales department staff are not only simply sales people but also engineers whose knowledge and experience would allow to find the best solution for each application. It is our philosophy. We never leave our customers without support. We approach them with great attention to appreciate their trust in us. Always there to help – we advise, train and make suggestions to our colleagues from the R&D department what needs arise in the market and what solutions the customers expect.

Małgorzata Szafarczyk, Head of Sales Department



### STANDARD FEATURES

- temperature range: +5°C above ambient temperature... +250°C
- other features like for drying ovens SL (see page 52)

### **ADVANTAGES OF SR HOT-AIR STERILIZERS**

- pre-set sterilization programs (including mask sterilization program)
- automatic door lock during the sterilization program
- automatically closed air-flap after starting the sterilization program
- 5 user programs and 3 pre-set programs

### AVAILABLE VERSIONS

- Smart
- Pass-through sterilizers

### SOFTWARE

 LabDesk for data download to a PC via LAN or Wi-Fi (optional for Smart version)



For surgical masks



For masks FFP2 and FFP3



For masks N95

# **Application**

- hot air sterilization
- disinfection of masks, documents, etc.




### SRW Smart hot-air sterilizers for mask disinfection


- Scientists from Stanford University recommend disinfecting N95 masks at 85°C for 20 minutes.
- The German government, in a published document, indicated recommendations for decontamination masks known as "surgical", FP2 and FFP3 masks by using hot (65-70°C) air for 30 minutes.

## SRW Smart hot-air sterilizers are an ideal solution for mask decontamination.

All of our units have a temperature display and a timer. They also have registration and full documentation of the process (duration and temperature).

Each unit also has a door lock – the lock is activated automatically when the program is started. This protects the whole process (all masks will complete the full disinfection cycle).





During one cycle you can decontaminate the following number of masks:

One cycle	Surgical masks	Masks FFP2, FFP3, N95	
SRW 115 Smart	150 pcs.	48 pcs.	
SRW 180 Smart	264 pcs.	80 pcs.	
SRW 240 Smart	300 pcs.	135 pcs.	



Sterilizers SR


		SR 53	SR 115	SR 240	SR 400	SR 750	SR 1000	
_		-	-				-41-	
Parameter								
air convection			natural (SRN) / forced (SRW)			forced (SRW)		
chamber capacity [I]		56	112	245	424	749	1005	
door type			solid/door with viewing window (option)					
temperature range			+5°C above ambient temperature+250°C					
temperature resolution [°C]				evei	ry 0,1			
controller		r	microprocessor PID, 4,3" (Smart) / 7" (Smart PRO) full colour touch screen					
interior			acid-proof stainless steel to DIN 1.4301					
housing	-			powder co	ated sheet			
	IG			stainless stee	el linen finish			
	A width	590	660	820	1020	1260	1260	
overall dims¹ [mm]	B height	710	850	1140	1430	1600	2000	
	C depth	620	710	770	770	880	880	
	D width	400	460	600	800	1040	1040	
internal dims [mm]	E height	390	540	800	1040	1200	1610	
	F depth	360	450	510	510	600	600	
30.3	-	25	25	25	25	-	-	
max shelf workload³ [kg]	PW <sup>2</sup> version	50	50	100	100	100	100	
max unit workload [kg]		40	60	90	120	140	300	
nominal power [W]		1700	2500	3100	4000	5500	5500	
weight <sup>s</sup> [kg]		48	65	114	162	260	307	
	SRN	0,4	0,4	0,6	-	-	-	
temperature fluctuation* at +105°C [+/- °C]	SRW	0,2	0,2	0,3	0,4	0,6	0,6	
	SRN	2,0	2,2	2,5	-	-	-	
temperature variation* at +105°C [+/- °C]	SRW	2,0	2,0	2,0	2,5	2,5	3,0	
over temperature protection			class 2.0 to DIN 12880 / class 3.1 (option)					
power supply**			230V 50-60Hz			400V 50-60Hz		
shelves fitted/max		2/5	2/7	3/10	3/14	5/16	6/22	
warranty			24 months					
manufacturer			POL-EKO-APARATURA					


all the above technical data refer to standard units (without optional accessories)

\* - fluctuation measured in centre of chamber; in space, variation (K) calculated for chamber as: K= +/- (T avg max - T avg min) / 2

\*\* - other power supplies on request

- 1 depth doesn't include 50 mm of power cable
- 2 reinforced shelf 3 on uniformly loaded surface





# Options and accessories (icon description see pages 76-82)




































	1		
-			

SRWP 240

Parameter			_			
air convection		forced	forced			
chamber capacity [l]		105	240			
door type		So	solid			
temperature range		+5°C above ambient	+5°C above ambient temperature+250°C			
temperature resolution [°C]		every 0,1				
controller		microprocessor PID, 4,3" (Smart) / 7'	microprocessor PID, 4,3" (Smart) /7" (Smart PRO) full colour touch screen			
interior		acid-proof stainles	acid-proof stainless steel to DIN 1.4301			
housing		powder co	pated sheet			
overall dims [mm]	A width	700	840			
	B height	910	1170			
	C depth	700	770			
internal dims [mm]	D width	460	600			
	E height	530	800			
	F depth	430	500			
max shelf workload [kg]		10	10			
PW version [kg]		50	100			
max unit worklad [kg]		60	90			
nominal power [W]		2500	3000			
weight [kg]		65	126			
over temperature protection		class 2.0 to DIN 12880 / class 3.1 (option)				
power supply**		230V 50-60Hz				
shelves fitted/max		2/7	3/10			
warranty		24 months				
manufacturer		POL-EKO-/	POL-EKO-APARATURA			

<sup>\*\* -</sup> other power supplies on request



# Options and accessories (icon description see pages 76-82)

























all the above technical data refer to standard units (without optional accessories)
\* - fluctuation measured in centre of chamber; in space, variation (K) calculated for chamber as: K = +/- (T avg max - T avg min) / 2